



# Configuration

## Projects

**Delete project**  
Note: All data must be removed from project first

**Rename project**

**Delete**

**Add/Edit**

**Edit settings**

**Members ?**

Name	Edit Data	Manage
Heather Arnaiz	✓	✓
Maria Edwards	<input type="checkbox"/>	<input type="checkbox"/> Remove

**Control permissions**

- Edit Data
  - Checked – user has vault role
  - Unchecked – read only or read export depending on vault role
- Manage
  - Checked – Manage users

**Select user from drop down**  
Note: Only members of a vault can be added to a project

- Select a user -  
- Select a user -  
Sylvia Ernst

## Members

**Your role**

- Click to view table of permissions for all roles

**Add a member**

- First name
- Last name
- Email
- Role

**Edit or delete member**

**Heather Arnaiz** Vault Administrator

**Add a member**

Name	Access	Date Added	Projects	Email	API Keys	Actions
Heather Arnaiz	vault admin					<input type="checkbox"/> <input type="checkbox"/>

First name: John  
Last name: Doe  
Email: john.doe@company.com  
Access role: Readonly

