

Vision

Filters



Vision

Plots

The screenshot displays the CDD.VISION interface with two plots and an export panel. The left plot is a scatterplot of log P vs. Molecular weight (g/mol). The right plot is a histogram of CNS MPO score. The export panel on the right allows for downloading plots as PDFs and includes options for axis labels, color, and size. Annotations highlight various features:

- Enlarge plot or table:** A magnifying glass icon over the plot area.
- Add selected molecules to collection:** A plus sign icon next to the 'Selected molecules 30 / 30' header.
- Export PDF of a chart:** The 'Export' button in the 'Export Plots' section.
- Drag to rearrange plots:** The 'Plots' section with a 'Drag the plots to reorder.' instruction and a visual representation of plot windows.

Export Plots
Download your plots as

Axis and Legend Labels
X: Molecular weight (g/mol)
Y: log P
Color: Single Point Screen: Avg Inhibiti...
Size: log D

Plots
Drag the plots to reorder.

Add a scatterplot or histogram to the display
Note: Maximum of 4 plots

Molecule	Molecular weight (g/mol)	log P	H-bond donors	H-bond acceptors	Lipinski violations	log D	log S
<input checked="" type="checkbox"/> CDDAT-0003209	154.125	0.23	1	5	0	-3.01	1

CC1=NC(=C(C=C1)C(=O)O)N